



MULTI-FUNCTION 4 IN 1 MACHINE

Descripción de productos

Real shot display



808 handle



Elight handle



ND YAG Laser handle



RF handle



Láser de diodo Elight/IPL Nd yag RF, 4 mangos diferentes en una sola máquina

1. Este 4 en 1 RF Elight IPL láser de diodo láser nd yag tiene accesorios de alta calidad, potencia fuerte * El mango del láser de diodo utiliza una barra láser dorada importada de Alemania, 50 millones de disparos.


*Lámpara del Reino Unido con mango Elight, 150 000-30 000 disparos

*Nd yag: Equipa 5 cabezales de tratamiento, ajustable 1064nm/532nm, fijo 1064nm/532nm, fijo 1320nm.


*RF: Equipa 6 cabezales de tratamiento, utilizados para la eliminación de arrugas, estiramiento de la piel, estiramiento de la piel.

Three wavelengths work together

12*24mm² {




12*35mm² {



**Suitable for hair removal on any part of the body,
Such as arms, legs, face, armpits, back, chest, etc.**





ELIGHT HANDLE

specification



- High-quality: Elight handpiece with super cooling system
- Optional Spot Size 15*50mm, 15*35mm
- Long shots time: 150000 times
- 7 filters: 420/510/560/585/640/690/755nm
Standard: 420/510/680nm



Before and after



Skin rejuvenation



Freckle removal



wrinkle reduction



Acne clearance



Hair removal



Pigmentation treatment

HANDLE-ADVANTAGE



specification

- Multiple choice
fixed 1064nm & 532nm, adjustable 1064nm & 532nm (adjust 1-7 mm) for tattoo removal, nevus removal etc.
- Multifunction
1320nm for skin rejuvenation
- Indicator light
Dot position is the indicator light position dot more accurate



Adjustable
532nm



Adjustable
1064nm



Black doll
1320nm



Fixed
532nm



Fixed
1064nm



Before and after



12mm eye-rf

28mm small tripolar rf head

36mm big tripolar rf head

14mm bipolar-rf

28mm bipolar-rf

38mm bipolar-rf

BEFORE & AFTER [®]



Hair removal



Skin rejuvenation



Freckle removal



Acne removal



Vascular



Lose weight

Advantage [®]



808 handle

Painless hair removal, Permanent hair removal, Suitable for all skin tones.

Elight handle

Elight for hair removal, skin rejuvenation, vascular therapy, pigmentation therapy, acne therapy.

ND YAG Laser handle

ND Yag laser for tattoo removal, eyeliner removal, eyebrow, embroidery removal, freckle removal.

RF handle

Skin tightening/lifting, wrinkle removal on body

**For more information about the introduction of the product and the details of use, contact us immediately
Online 24 hours!!**