

New upgrade

# 3IN1 HIFU

Macro & Micro High Intensity Focused Ultrasound System



SAFE



EFFECTIVE



MULTI-DEPTH



MMFU (Macro & Micro Focused Ultrasound) stands for "Macro & Micro High Energy Focused Ultrasound System" non-surgical treatment of face lifting, firming and body contouring system.



\*A high repetition rate of short pulses are delivered deep into the dermis, achieving high average power and therapeutically effective heat build-up, with virtually no pain. \*The latest machine has total 11 probes, which are fully carved and anti-aging.

\*Non surgical treatment of facial lifting, tightening and body contour shaping system. The world's latest high energy focused ultrasound technology.


# 4K HD touch screen

Experience the World in 4K HD  
Every Detail, Every Moment, Every Color



## Color touch screenskin rejuvenation

1. 15.6 inch big touch screen, humanized designed
2. Male/femal, body parts, 6 skin types for option.
3. English+your native language.
4. Free OEM program service(logo/system...)









## The Vagina Tightening HIFU System

- Vaginal fast tightening elasticity enhancement
- Vaginal tissue collagen rebirth
- Eliminate dry pinch and improve lubrication
- Alleviate menopausal symptoms
- Non-invasive, painless, no recovery
- Safe and reliable
- Using world-class high energy (focusing ultrasound body technology)
- Better rejuvenation than the laser



## For whole body

1.	For fine lines		1.5mm 7MHz
2.	For cuticle		2.0mm 5.5MHz
3.	For dermis		3.0mm 7MHz
4.	For fascial layer		4.5mm 4MHz
5.	For double chins		6.0mm 2MHz
6.	For arms		9.0mm 2MHz
7.	For the hips and waistline		13mm 2MHz
8.	For fascia layer		1.5mm 7MHz
9.	For the arm bye meat		3.0mm 7MHz
10.	For the hips and waistline		4.5mm 4MHz
11.	Private treatment head		3.0mm 4.5mm



## Comparison before and after use



1. Remove wrinkles on around forehead, eyes, mouth, etc.
2. Lifting and tightening both cheeks skin.
3. Improving skin elasticity and shaping contour.
4. Improving jaw line, reducing "marionette lines".
5. Tightening the skin tissue on forehead, lifting the eyebrows lines.
6. Improving skin complexion, making the skin delicate and bright.
7. Match with injecting beauty like Hyaluronicacid, collagen, to resolve more aging problem.
8. Removing neck wrinkles, protecting neck aging.



---

## Product Specifications

<b>Device Name</b>	3IN1 HIFU
<b>Technology</b>	HIFU
<b>Input Voltage</b>	AC100V - 240V
<b>Voltage</b>	2-7MHz
<b>Maximum Power</b>	200-400VA
<b>Fuse</b>	5A
<b>Number Of Handles</b>	3 Handles
<b>Treatment Area</b>	Face Body
<b>Cartridges</b>	N1.5mm, S2.0mm, N3.0mm, N4.5mm, W6.0mm, W9.0mm, W13.0mm, P1.5, P3.0, P4.5 DS7-3.0, DS7-4.5
<b>Shots</b>	12000 Shots / Cartridge
<b>Total energy</b>	10 to 120mi/pin adjustable
<b>Packing Weight</b>	12.3kg
<b>Packing Size</b>	556*570*235mm

\*A high repetition rate of short pulses are delivered deep into the dermis, achieving high average power and therapeutically effective heat build-up, with virtually no pain. \*The latest machine has total 11 probes,



which are fully carved and anti-aging.

\*Non surgical treatment of facial lifting, tightening and body contour shaping system. The world's latest high energy focused ultrasound technology.